

Swing™

BREASTPUMP



The only single electric, personal use pump with the breakthrough 2-Phase Expression™ Technology, designed to get more milk in less time*. First researched and developed for hospital breastpumps, 2-Phase Expression is an advanced pumping pattern that mimics the baby's natural pattern.

HOW DO BABIES NURSE?

*Studies have shown**, that babies nurse in two phases:*

- Fast and light sucking to start the milk flowing.
- When the milk starts flowing the baby switches to slower, deeper (stronger) sucks, to get milk as quickly as possible.

2-Phase Expression technology has been designed to mimic a baby's natural nursing rhythm.

- Stimulation phase (*phase 1*) fast and light to start the milk flowing.
- Expression phase (*phase 2*) slower and with more vacuum to express more milk gently and efficiently.

*Swing,
only single electric
Breastpump
with
2-Phase Expression*



Swing Breastpump
#67050

*When pumping at Maximum Comfort Vacuum

**Woolridge MW, The 'Anatomy' of Instant Sucking *Midwifery*, (1986) 2, 164-171



Swing™

BREASTPUMP

For the mom on the go!

- One touch let-down button designed for faster milk flow
- Convenient – easily fits into a purse or briefcase
- Efficient, adjustable speed and vacuum
- Comfortable SoftFit™ breastshield
- Mobile, use on a tabletop, belt clip or shoulder/neck strap



Swing Breastpump
#67050



Item	Article #	UPC Code	Unit Count	Unit Dimensions HxWxD - Inches	Unit Wt. Lbs.	Case Pack/ Carton	Shipping Carton Dimensions HxWxD - Inches	Shipping Carton Weight - Lbs.
Swing Single Electric	67050	0 20451 67050 8	1	9 ^{5/8} x 10 ^{7/8} x 4 ^{11/16}	2 lbs. 2oz.	2	10 ^{1/8} x 10 x 11 ^{5/8}	5 lbs.



Medela, Inc., P.O. Box 660, McHenry, IL 60051-0660 USA
Phone: 888-MEDELA-8 (888-633-3528) or 815-363-1166
Fax: 815-363-1246
Email: customer.service@medela.com www.medela.com

Swing and SoftFit are trademarks and 2-Phase Expression and Medela are registered trademarks of Medela, Inc. ©2006 Medela, Inc.
1907732 A 0406 © 2006 Medela, Inc. Printed in USA

To order or for more information, call your Medela Customer Service team today at 888-MEDELA-8 (888-633-3528) or visit www.medela.com